

Human Acute Kidney Injury Antibody Array - Membrane (20 targets) ab169806

1 References [画像数 1](#)

製品の概要

製品名	Human Acute Kidney Injury Antibody Array - Membrane (20 targets)
サンプルの種類	Cell culture supernatant, Saliva, Milk, Urine, Serum, Plasma, Cell culture extracts, Other biological fluids, Whole Blood, Tissue Extracts, Cell Lysate, Cell culture media
アッセイタイプ	Semi-quantitative
種交差性	交差種: Human
製品の概要	<p>For the simultaneous detection of 20 Human Cytokine levels in 4 or 8 samples. Suitable for all sample types.</p> <p>Targets:</p> <p>KIM-1, Albumin (ALB), Osteopontin (Osteopontin), Trefoil Factor 3 (TFF3), B2M, Clusterin (CLU), CXCL16, Osteoactivin, Liver-type fatty acid binding protein (L-FABP), MCP-1, sTNFR1, Calbindin (CALB1), Interferon gamma-induced protein 10 (IP-10), Cystatin C, HGF, MIF, NGAL, TIMP-1, VCAM-1, VEGF-A</p> <p><u>Learn more about membrane antibody arrays</u></p>

特記事項	<p>Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances.</p> <p>It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.</p>
------	--

アプリケーション	適用あり: Multiplex Protein Detection
----------	--

製品の特性

保存方法	Store at -20°C. Please refer to protocols.
------	--

内容	1 x 4 Membranes	1 x 8 Membranes
1,000X HRP-Conjugated Streptavidin	1 x 50µl	1 x 50µl
1X Blocking Buffer	1 x 25ml	2 x 25ml


内容	1 x 4 Membranes	1 x 8 Membranes
2000X Biotin-Conjugated Anti-Cytokines	2 vials	4 vials
20X Wash Buffer I	1 x 10ml	1 x 20ml
20X Wash Buffer II	1 x 10ml	1 x 20ml
2X Cell Lysis Buffer	1 x 10ml	1 x 16ml
8-Well Incubation Tray (with Lid)	1 unit	1 unit
Detection Buffer C	1 x 1.5ml	1 x 2.5ml
Detection Buffer D	1 x 1.5ml	1 x 2.5ml
Human Acute Kidney Injury Antibody Array Membranes	4 units	8 units

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab169806の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Multiplex Protein Detection		Use at an assay dependent concentration.

画像



Typical results obtained with Abcam Membrane Antibody Arrays.
Image above presents typical images obtained with Abcam Human Cytokine Membrane Antibody Array. These membranes were probed with conditioned media from two different cell lines.
Membranes were exposed to film at RT for 1 min.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors