

# Human Cytokine Antibody Array - Membrane (60 targets) ab169818

**1 References**    [画像数 1](#)

### 製品の概要

製品名	Human Cytokine Antibody Array - Membrane (60 targets)
サンプルの種類	Cell culture supernatant, Saliva, Milk, Urine, Serum, Plasma, Cell culture extracts, Other biological fluids, Whole Blood, Tissue Extracts, Cell Lysate, Cell culture media
アッセイタイプ	Semi-quantitative
種交差性	交差種: Human
製品の概要	ab169818 is for simultaneous detection of 60 Human Cytokines in 4 or 8 samples. Suitable for all liquid sample types.

**Targets:** Adiponectin, AgRP, Angiopoietin 2, Amphiregulin, Axl, bFGF, b-NGF, BTC, CCL28, CTACK, Dtk, EGFR, ENA-78, Fas, TNFRSF6, FGF-4, FGF-9, GCSF, GTR-Ligand, GTR, GRO, GRO alpha, HCC-4, HGF, ICAM-1, ICAM-3, IGFBP-3, IGFBP-6, IGF-1 sR, IL-1 R4/ST2, IL-1 RI, IL-11, IL-12 p40, IL-12 p70, IL-17, IL-2 ra, IL-6 R, IL-8, I-TAC, Lymphotoxin, MIF, MIP-1 alpha, MIP-1 beta, MIP-3 beta, MSP-alpha chain, NT-4, Osteoprotegerin, Oncostatin M, PLGF, sgp130, sTNFRII, sTNFRI, TECK, TIMP-1, TIMP-2, Thrombopoietin, TRAIL R3, TRAIL R4, uPAR, VEGF-A, VEGF-D

### [Learn more about membrane antibody arrays](#)

特記事項	Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances. It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.
------	---

アプリケーション	<b>適用あり:</b> Multiplex Protein Detection
----------	--

### 製品の特性

保存方法	Store at -20°C. Please refer to protocols.
------	--

内容	1 x 4 Membranes	1 x 8 Membranes	1 x 2 Membranes
1,000X HRP-Conjugated Streptavidin	1 x 50µl	1 x 50µl	1 x 50µl


内容	1 x 4 Membranes	1 x 8 Membranes	1 x 2 Membranes
1X Blocking Buffer	1 x 25ml	2 x 25ml	1 x 25ml
20X Wash Buffer I	1 x 10ml	1 x 20ml	1 x 10ml
20X Wash Buffer II	1 x 10ml	1 x 20ml	1 x 10ml
2X Cell Lysis Buffer	1 x 10ml	1 x 16ml	1 x 10ml
8-Well Incubation Tray (with Lid)	1 unit	1 unit	1 unit
Biotinylated Antibody Cocktail (C7)	2 vials	4 vials	1 vial
Detection Buffer C	1 x 1.5ml	1 x 2.5ml	1 x 1.5ml
Detection Buffer D	1 x 1.5ml	1 x 2.5ml	1 x 1.5ml
Human Cytokine Antibody Array Membranes (C7)	4 units	8 units	2 units

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab169818の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
<b>Multiplex Protein Detection</b>		Use at an assay dependent concentration.

## 画像



Typical images obtained with Abcam Human Cytokine Antibody Array. The membranes were probed with conditioned media from two different cell lines or a mock buffer control. Membranes were exposed to film at room temperature for 1 minute.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors