# abcam

### Product datasheet

## Alexa Fluor® 488 Anti-HDAC2 antibody [Y461] ab196471



ועלטעבעו RabMAb

#### 画像数 4

#### 製品の概要

免疫原

特記事項

製品名 Alexa Fluor® 488 Anti-HDAC2 antibody [Y461]

製品の詳細 Alexa Fluor® 488 Rabbit monoclonal [Y461] to HDAC2

由来種 Rabbit

Alexa Fluor® 488. Ex: 495nm, Em: 519nm 標識

アプリケーション 適用あり: ICC/IF, Flow Cyt (Intra)

種交差性 交差種: Human

交差が予測される動物種: Mouse, Rat 4

ポジティブ・コントロール ICC/IF: HeLa and wildtype HAP1 cells. Flow Cyt (intra): HeLa cells.

> Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb® patents.

Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.

Alexa Fluor® is a registered trademark of Molecular Probes, Inc, a Thermo Fisher Scientific Company. The Alexa Fluor® dye included in this product is provided under an intellectual property license from Life Technologies Corporation. As this product contains the Alexa Fluor® dye, the purchase of this product conveys to the buyer the non-transferable right to use the purchased product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). As this product contains the Alexa Fluor® dye the sale of this product is expressly conditioned on the buyer not using the product or its components, or any materials made using the product or its components, in any activity to generate revenue, which may include, but is not limited to use of the product or its components: in manufacturing; (ii) to provide a service, information, or data in return for payment (iii) for therapeutic, diagnostic or prophylactic purposes; or (iv) for resale, regardless of whether they are sold for use in research. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or

outlicensing@thermofisher.com.

#### 製品の特性

製品の状態

保存方法 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C.

Avoid freeze / thaw cycle. Stable for 12 months at -20°C. Store In the Dark.

**バッファー** pH: 7.40

Preservative: 0.02% Sodium azide

Constituents: PBS, 30% Glycerol (glycerin, glycerine), 1% BSA

精製度 Protein A purified

**ポリ/モノ** モノクローナル

**クローン名** Y461 アイソタイプ IgG

#### アプリケーション

**The Abpromise guarantee** <u>Abpromise保証は、</u>次のテスト済みアプリケーションにおけるab196471の使用に適用されます アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/100 - 1/500. This product gave a positive signal in both 100% methanol (5 min) and 4% formaldehyde (10min) fixed HeLa cells.
Flow Cyt (Intra)		1/500.  ab199091 - Rabbit monoclonal IgG (Alexa Fluor® 488), is suitable for use as an isotype control with this antibody.


#### ターゲット情報

機能	Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones	
	(H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an	
	important role in transcriptional regulation, cell cycle progression and developmental events.	
	Histone deacetylases act via the formation of large multiprotein complexes.	
	Forms transcriptional repressor complexes by associating with MAD, SIN3, YY1 and N-COR.	
	Interacts in the late S-phase of DNA-replication with DNMT1 in the other transcriptional repressor	
	complex composed of DNMT1, DMAP1, PCNA, CAF1. Deacetylates TSHZ3 and regulates its	
	transcriptional repressor activity.	
組織特異性	Widely expressed; lower levels in brain and lung.	
配列類似性	Belongs to the histone deacetylase family. HD type 1 subfamily.	
翻訳後修飾	S-nitrosylated by GAPDH. In neurons, S-Nitrosylation at Cys-262 and Cys-274 does not affect the	

S-nitrosylated by GAPDH. In neurons, S-Nitrosylation at Cys-262 and Cys-274 does not affect the enzyme activity but abolishes chromatin-binding, leading to increases acetylation of histones and activate genes that are associated with neuronal development. In embryonic cortical neurons, S-Nitrosylation regulates dendritic growth and branching.

細胞内局在 Nucleus.


#### 画像



Immunocytochemistry/ Immunofluorescence - Alexa Fluor® 488 Anti-HDAC2 antibody [Y461] (ab196471)

ab196471 staining HDAC2 in wild-type HAP1 cells (top panel) and HDAC2 knockout HAP1 cells (bottom panel). The cells were fixed with 100% methanol (5min), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with ab196471 at 1/500 dilution (shown in green) and ab195889 at 1/250 dilution (shown in pseudocolour red) overnight at +4°C. Nuclear DNA was labelled in blue with DAPI.


Image was taken with a confocal microscope (Leica-Microsystems, TCS SP8).



Flow Cytometry (Intracellular) - Alexa Fluor® 488 Anti-HDAC2 antibody [Y461] (ab196471)

Overlay histogram showing HeLa cells stained with ab196471 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab196471, 1/500 dilution) for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) Alexa Fluor® 488 used at the same concentration and conditions as the primary antibody. Unlabelled sample (blue line) was also used as a control.

Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter.




Immunocytochemistry/ Immunofluorescence - Alexa Fluor® 488 Anti-HDAC2 antibody [Y461] (ab196471)

ab196471 staining HDAC2 in HeLa cells. The cells were fixed with 100% methanol (5 min), permeabilized in 0.1% Triton X-100 for 5 minutes and then blocked in 1% BSA/10% normal goat serum/0.3M glycine in 0.1%PBS-Tween for 1h. The cells were then incubated with ab196471 at a working dilution of 1 in 100 (shown in green) and **ab195889**, Mouse monoclonal [DM1A] to alpha Tubulin (Alexa Fluor® 594, shown in red) at 2µg/ml overnight at +4°C. Nuclear DNA was labelled in blue with DAPI.

This product gave a positive signal in 4% formaldehyde (10 min) fixed HeLa cells under the same testing conditions.

Image was taken with a confocal microscope (Leica-Microsystems, TCS SP8).



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <a href="https://www.abcam.co.jp/abpromise">https://www.abcam.co.jp/abpromise</a> or contact our technical team.

#### Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors