abcam

Product datasheet

FITC Anti-Hemoglobin antibody ab19361

製品の概要

製品名 FITC Anti-Hemoglobin antibody

製品の詳細 FITC Goat polyclonal to Hemoglobin

由来種 Goat

標識 FITC. Ex: 493nm, Em: 528nm

アプリケーション 適用あり: Immunoelectrophoresis, ELISA, ICC

種交差性 交差種: Human

交差が予測される動物種: Chimpanzee, Rhesus monkey, Gorilla 🔷

免疫原 Full length protein corresponding to Human Hemoglobin. Goats were immunized with purified

human hemoglobin A1 of adult red blood cell origin. T

特記事項 F/P ratio is 3 to 7.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets

your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be

found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Store at +4°C.

バッファー pH: 6.8

Preservative: 0.1% Sodium azide Constituents: PBS, 0.2% BSA

精製度 Immunogen affinity purified

特記事項(精製) The antibody was isolated by affinity chromatography using antigen coupled to agarose beads.

Antibody concentration was determined by extinction coefficient prior to conjugation: absorbance

at 280 nm of 1.4 equals 1.0 mg of lgG.

ポリ/モノ ポリクローナル

アイソタイプ lgG

1

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab19361の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Immunoelectrophoresis		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
ICC		Use at an assay dependent concentration.

ターゲット情報

機能

組織特異性

関連疾患

Involved in oxygen transport from the lung to the various peripheral tissues.

Red blood cells.

Defects in HBA1/HBA2 may be a cause of Heinz body anemias (HEIBAN) [MIM:140700]. This is a form of non-spherocytic hemolytic anemia of Dacie type 1. After splenectomy, which has little benefit, basophilic inclusions called Heinz bodies are demonstrable in the erythrocytes. Before splenectomy, diffuse or punctate basophilia may be evident. Most of these cases are probably instances of hemoglobinopathy. The hemoglobin demonstrates heat lability. Heinz bodies are observed also with the Ivemark syndrome (asplenia with cardiovascular anomalies) and with glutathione peroxidase deficiency.

Defects in HBA1/HBA2 are the cause of alpha-thalassemia (A-THAL) [MIM:604131]. The thalassemias are the most common monogenic diseases and occur mostly in Mediterranean and Southeast Asian populations. The hallmark of alpha-thalassemia is an imbalance in globin-chain production in the adult HbA molecule. The level of alpha chain production can range from none to very nearly normal levels. Deletion of both copies of each of the two alpha-globin genes causes alpha(0)-thalassemia, also known as homozygous alpha thalassemia. Due to the complete absence of alpha chains, the predominant fetal hemoglobin is a tetramer of gamma-chains (Bart hemoglobin) that has essentially no oxygen carrying capacity. This causes oxygen starvation in the fetal tissues leading to prenatal lethality or early neonatal death. The loss of three alpha genes results in high levels of a tetramer of four beta chains (hemoglobin H), causing a severe and lifethreatening anemia known as hemoglobin H disease. Untreated, most patients die in childhood or early adolescence. The loss of two alpha genes results in mild alpha-thalassemia, also known as heterozygous alpha-thalassemia. Affected individuals have small red cells and a mild anemia (microcytosis). If three of the four alpha-globin genes are functional, individuals are completely asymptomatic. Some rare forms of alpha-thalassemia are due to point mutations (non-deletional alpha-thalassemia). The thalassemic phenotype is due to unstable globin alpha chains that are rapidly catabolized prior to formation of the alpha-beta heterotetramers.

Note=Alpha(0)-thalassemia is associated with non-immune hydrops fetalis, a generalized edema of the fetus with fluid accumulation in the body cavities due to non-immune causes. Non-immune hydrops fetalis is not a diagnosis in itself but a symptom, a feature of many genetic disorders, and the end-stage of a wide variety of disorders.

配列類似性

翻訳後修飾

Belongs to the globin family.

The initiator Met is not cleaved in variant Thionville and is acetylated.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors