# abcam

# Product datasheet

# Anti-Glycine antibody ab9442

## ★★★★★ 1 Abreviews 9 References

#### 製品の概要

製品名 Anti-Glycine antibody

製品の詳細 Rabbit polyclonal to Glycine

由来種 Rabbit

特異性 The antibody is calibrated against a spectrum of antigens to assure hapten selectivity. No

measurable cross-reactivity (<1:1000) was detected against glycine in peptides or proteins. Fixed tissue cross-reactivity was tested with known targets at the recommended dilution. No measurable glutaraldehyde-fixed tissue cross-reactivity (<1:1000) was detected against L-alanine, gamma-aminobutyrate, agmatine, guanidine, D/L-arginine, L-citrulline, L-cysteine, D/L-glutamate, D/L-glutamine, glutathione, L-lysine, L-ornithine, L-serine, taurine, L-threonine, L-

tryptophan, L-tyrosine.

アプリケーション 適用あり: IHC-FoFr, IHC-Fr, Immunomicroscopy

種交差性 交差種: Species independent

免疫原 Chemical/ Small Molecule corresponding to Glycine (Glutaraldehyde).

特記事項 The product is optimized for HPI/EHPI with gold or fluorescence detection using etched plastic

sections. Filter diluted reagents with 0.2 mm syringe filters before use on EM grids. Enzymelinked visualizations can be used but will compress the signal dynamic range and are less

sensitive.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years.

Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets

your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be

found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze /

thaw cycle.

精製度 IgG fraction ポリクローナル

1

#### アプリケーション

**The Abpromise guarantee** <u>Abpromise保証は、</u>次のテスト済みアプリケーションにおけるab9442の使用に適用されます アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-FoFr		Use at an assay dependent concentration. Use with frozen or vibratome sections is possible but will not yield optimal images as IgGs penetrate aldehyde cross-linked tissue poorly and most amino acids are present at such high levels that prozone effects occur. Use in whole mounts is not recommended for similar reasons.
IHC-Fr	**** (1)	Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.

#### ターゲット情報

#### 関連性

Defects in GLDC are a cause of nonketotic hyperglycinemia (NKH), also known as glycine encephalopathy (GCE). NKH is an autosomal recessive disease characterized by accumulation of a large amount of glycine in body fluid and by severe neurological symptoms. The degredation of glycine is catalised by the glycine cleavage system. The P protein binds the alpha-amino group of glycine through its pyridoxal phosphate cofactor; carbondioxide is released and the remaining methylamine moiety is then transferred to the lipoamide cofactor of the H protein. The glycine cleavage system is composed of four proteins: P, T, L and H.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <a href="https://www.abcam.co.jp/abpromise">https://www.abcam.co.jp/abpromise</a> or contact our technical team.

#### Terms and conditions

•	Guarantee only valid for products bought direct from Abcam or one of our authorized distributors				
		3			