

# Anti-Marek's Disease antibody ab33718

### 製品の概要

製品名	Anti-Marek's Disease antibody
製品の詳細	Chicken polyclonal to Marek's Disease
由来種	Chicken
特異性	This antibody is monospecific for Marek s Disease.
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Gallid alphaherpesvirus 2
免疫原	Tissue, cells or virus corresponding to Chicken Marek's Disease.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Constituents: Whole serum, 40% Glycerol
精製度	Whole antiserum
特記事項(精製)	Sterile filtered serum produced in SPF (Specific Pathogen Free) chicken embryos.
ポリモノ	ポリクローナル
アイソタイプ	IgY

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab33718の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.

## ターゲット情報

---

**関連性** Marek's Disease Virus, (MDV), is an oncogenic avian double stranded DNA herpesvirus. It is highly contagious and results in aggressive T-cell lymphomas. Infection of MDV occurs via the respiratory tract, where its uptake by macrophages and dendritic cells allows it to be transported in the bloodstream to visceral organs and peripheral nerves, causing death in up to 80 percent of birds. It occurs in two main forms, depending on which parts of the body are affected by the tumours. In the nervous form the nerves, particularly the sciatic nerves, are affected. In the visceral form, greyish-white tumours are found in the ovaries, liver, spleen, kidney, heart and other organs.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors