abcam

Product datasheet

Human ABCB7 peptide ab66358

製品の詳細

製品名 Human ABCB7 peptide

精製度 > 70 % HPLC.

70 - 90% by HPLC

Animal free No

由来 Synthetic

生物種 Human

特性

Our Abpromise guarantee covers the use of ab66358 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking

製品の状態 Liquid

備考 - First try to dissolve a small amount of peptide in either water or buffer. The more charged

residues on a peptide, the more soluble it is in aqueous solutions.

- If the peptide doesn't dissolve try an organic solvent e.g. DMSO, then dilute using water or

buffer.

- Consider that any solvent used must be compatible with your assay. If a peptide does not

dissolve and you need to recover it, lyophilise to remove the solvent.

- Gentle warming and sonication can effectively aid peptide solubilisation. If the solution is

cloudy or has gelled the peptide may be in suspension rather than solubilised.

- Peptides containing cysteine are easily oxidised, so should be prepared in solution just prior

to use.

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

Information available upon request.

関連情報

1

機能 Could be involved in the transport of heme from the mitochondria to the cytosol. Plays a central

role in the maturation of cytosolic iron-sulfur (Fe/S) cluster-containing proteins.

関連疾患 Defects in ABCB7 are the cause of X-linked sideroblastic anemia with ataxia (ASAT)

[MIM:301310]. ASAT is a recessive disorder characterized by an infantile to early childhood onset

of nonprogressive cerebellar ataxia and mild anemia with hypochromia and microcytosis.

配列類似性 Belongs to the ABC transporter superfamily. ABCB family. Heavy Metal importer (TC 3.A.1.210)

subfamily.

Contains 1 ABC transmembrane type-1 domain.

Contains 1 ABC transporter domain.

細胞内局在 Mitochondrion inner membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors