abcam

Product datasheet

Recombinant Human SPR protein ab86445

画像数1

製品の詳細

製品名 Recombinant Human SPR protein

精製度 > 95 % SDS-PAGE.

ab86445 is purified using conventional chromatography techniques.

発現系 Escherichia coli

タンパク質長 Full length protein

Animal free No.

由来 Recombinant

生物種 Human

配列 MGSSHHHHHH SSGLVPRGSH MEGGLGRAVC

LLTGASRGFG RTLAPLLASL LSPGSVLVLS
ARNDEALRQL EAELGAERSG LRVVRVPADL
GAEAGLQQLL GALRELPRPK GLQRLLLINN
AGSLGDVSKG FVDLSDSTQV NNYWALNLTS
MLCLTSSVLK AFPDSPGLNR TVVNISSLCA
LQPFKGWALY CAGKAARDML FQVLALEEPN
VRVLNYAPGP LDTDMQQLAR ETSVDPDMRK
GLQELKAKGK LVDCKVSAQK LLSLLEKDEF

KSGAHVDFYD K

特性

Our **Abpromise guarantee** covers the use of **ab86445** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

備考 Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method).

前処理および保存

保存方法および安定性 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

1

関連情報

機能 Catalyzes the final one or two reductions in tetra-hydrobiopterin biosynthesis to form 5,6,7,8-

tetrahydrobiopterin.

関連疾患 Defects in SPR are the cause of dystonia DOPA-responsive due to sepiapterin reductase

deficiency (DRDSPRD) [MIM:612716]. In the majority of cases, patients manifest progressive psychomotor retardation, dystonia and spasticity. Cognitive anomalies are also often present. The disease is due to severe dopamine and serotonin deficiencies in the central nervous system caused by a defect in BH4 synthesis. Dystonia is defined by the presence of sustained involuntary

muscle contractions, often leading to abnormal postures.

配列類似性 Belongs to the sepiapterin reductase family.

翻訳後修飾 In vitro phosphorylation of Ser-213 by CaMK2 does not change kinetic parameters.

細胞内局在 Cytoplasm.

画像



15% SDS-PAGE showing ab86445 at approximately 30.2kDa (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors