abcam

Product datasheet

PAP pen ab2601

★★★★★ 4 Abreviews 19 References 画像数 1

製品の概要

製品名

特記事項

PAP pen

The PAP pen is a special marking pen that provides a thin film-like green-tinged hydrophobic barrier when a circle is drawn gently around a specimen on a slide. Gentle application of the pen tip while drawing is advised.

This water repellent barrier keeps staining reagents localized on the tissue sections, prevents mixing of reagents when differentially staining two sections on the same slide, and allows use of less reagents per section.

Insoluble in alcohol and acetone, PAP Pen markings can be optionally removed with xylene after the staining procedure is complete.

PAP Pen is good for approximately 500 uses.

The formulation of PAP pen has been optimised to increase water repellence and tolerance against alcohol and surface active agent, and is now a lighter colour.

製品の特性

製品の状態

Liquid

保存方法

Store at room temperature.

バッファー

pH: 5.8

Constituents: 35% Toluene, 15% Solvent naphtha, 15% Xylene

画像



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors