

Phosphatase Inhibitor Cocktail II ab201113

★★★★★ [1 Abreviews](#) [2 References](#)

製品の概要

製品名	Phosphatase Inhibitor Cocktail II
製品の概要	<p>Abcam's phosphatase inhibitor cocktails inhibit a number of phosphatases to preserve the phosphorylation state of the proteins in your sample.</p> <p>Target phosphatases - alkaline phosphatases, acid phosphatases, serine/threonine phosphatases, protein tyrosine phosphatases</p> <p>Inhibitors included - imidazole, sodium fluoride, sodium molybdate, sodium orthovanadate, sodium tartrate dihydrate</p> <p>Form: Supplied in water (1 ml)</p> <p>Preparation: The recommended starting dilution is 1:100</p> <p>Visit our inhibitor cocktails page for more information on our complete range of convenient and easy to use protease and phosphatase inhibitor cocktails.</p>
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called K275 EZBlock™ Phosphatase Inhibitor Cocktail II. K275-1 is the same size as the 1 ml size of ab201113.</p> <p>Individual components and concentrations:</p> <p>Imidazole: alkaline phosphatase inhibitor (200 mM)</p> <p>Sodium fluoride: acid phosphatase and serine/threonine phosphatase inhibitor (100 mM)</p> <p>Sodium molybdate: acid phosphatase inhibitor (115 mM)</p> <p>Sodium orthovanadate: protein tyrosine phosphatase and alkaline phosphatase inhibitor (100 mM)</p> <p>Sodium tartrate dihydrate: acid phosphatase inhibitor (400 mM)</p>

製品の特性

保存方法 Store at -20°C. Please refer to protocols.

内容	1 ml
Phosphatase Inhibitor Cocktail II	1 x 1ml

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors