

# Trypsin Antigen Retrieval Solution ab970

★★★★★ **1 Abreviews** **24 References**

### 製品の概要

製品名	Trypsin Antigen Retrieval Solution
アプリケーション	<b>適用あり:</b> IHC-P
特記事項	<p>Trypsin-based antigen retrieval solution for Protease-Induced Epitope Retrieval (PIER) (ab970). Contains colored pH indicator for easy pH confirmation.</p> <p>Protocol Notes: 1) Deparaffinize tissues and hydrate to water. If necessary perform a hydrogen peroxide block, wash in water, and rinse in PBS or TBS. The tissue section is ready for trypsin digestion. 2) Combine 1 part trypsin concentrate with 1 part buffer (working solution stable for 5 working days at 2-8°C). Incubate section for 5-10 minutes at 37°C or 15 minutes at room temperature. Then, rinse well in PBS or TBS buffer and proceed with immunostaining procedure. This product utilises a color-coded pH indicator (rose) for easy identification. The end-user can visually inspect the solution and see that the mixture or buffer is at the proper pH. If the mixed solution is purple, the pH is too high for optimum digestion. If the buffer or trypsin solution turns orange or yellow, the pH is too low for optimum digestion.</p> <p>Trypsin is a commonly used digestive enzyme. In formalin-fixed, paraffin-embedded tissues, certain antibody protocols required enzyme pretreatment for proper immunohistochemical staining.</p> <p><b>Other products for IHC</b></p> <p>Other enzymes for antigen retrieval include: <b><u><a href="#">Proteinase K Antigen Retrieval Solution ab64220</a></u></b> and <b><u><a href="#">Pepsin Solution ab64201</a></u></b>.</p> <p>Find more kits and reagents for antigen retrieval, blocking, signal amplification, visualization, counterstaining, and mounting in the <b><u><a href="#">IHC kits and reagents guide</a></u></b>.</p>

### 製品の特性

製品の状態	Liquid
保存方法	Store at +4°C.
バッファー	Preservative: 0.05% Sodium azide Constituents: PBS, Carrier protein
IHC 試薬 備考	Trypsin is a commonly used digestive enzyme. In formalin-fixed, paraffin-embedded tissues, certain antibody protocols required enzyme pretreatment for proper immunohistochemical staining.

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab970の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent dilution. Ratio for 25ml concentrated trypsin and 25 ml buffer is 1:1 - 1:3.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors