

# Human Fetuin-A Antibody Pair - BSA and Azide free ab270374

リコンビナント **RabMAb**

画像数 4

## 製品の概要

|            |   |
|------------|---|
| 製品名        | Human Fetuin-A Antibody Pair - BSA and Azide free   |
| アッセイタイプ    | ELISA set   |
| 検出範囲       | 203.125 pg/ml - 13000 pg/ml   |
| 種交差性       | 交差種: Human  |
| 製品の概要      | <p>The Antibody Pair can be used to quantify Human Fetuin-A. BSA and Azide free antibody pairs include unconjugated capture and detector antibodies suitable for sandwich ELISAs. The antibodies are provided at an approximate concentration of 1 mg/ml as measured by the protein A280 method. The recommended antibody orientation is based on internal optimization for ELISA-based assays. Antibody orientation is assay dependent and needs to be optimized for each assay type. Both capture and detector antibodies are rabbit monoclonal antibodies delivering consistent, specific, and sensitive results.</p> <p>For additional information on the performance of the antibody pair, see the equivalent SimpleStep ELISA® Kit (<a href="#">ab269372</a>), which uses the same antibodies. However, due to differences in their formulation, this antibody pair cannot be used with the consumables provided with our SimpleStep ELISA Kits. Please note that the range provided for the pairs is only an estimation based on the performance of the related product using the same antibody pair. Performance of the antibody pair will depend on the specific characteristics of your assay. We guarantee the product works in sandwich ELISA, but we do not guarantee the sensitivity or dynamic range of the antibody pair in your assay.</p> <p>Download SDS <a href="#">here</a>.</p> |
| アプリケーション   | 適用あり: Sandwich ELISA  |
| 試験プラットフォーム | Reagents  |

## 製品の特性

|          |   |
|----------|---|
| 保存方法     | Store at +4°C. Please refer to protocols. |
| キャリア・フリー | はい  |

| 内容  | 10 x 96 tests |
|---|---------------|
| Human Fetuin-A Capture Antibody (unconjugated)  | 1 x 100µg     |
| Human Fetuin-A Detector Antibody (unconjugated) | 1 x 100µg     |


|       |  |
|-------|--|
| 機能    | Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.                               |
| 組織特異性 | Synthesized in liver and selectively concentrated in bone matrix. Secreted in plasma. It is also found in dentin in much higher quantities than other plasma proteins. |
| 配列類似性 | Belongs to the fetuin family.<br>Contains 2 cystatin fetuin-A-type domains.  |
| 翻訳後修飾 | Phosphorylation sites are present in the extracellular medium.   |
| 細胞内局在 | Secreted.  |

## アプリケーション


**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab270374の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション              | Abreviews | 特記事項                                     |
|-----------------------|-----------|--|
| <b>Sandwich ELISA</b> |           | Use at an assay dependent concentration. |

## 画像




Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab269372](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.



Sandwich ELISA - Human Fetuin-A Antibody Pair -  
BSA and Azide free (ab270374)

Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab269372](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.



Sandwich ELISA - Human Fetuin-A Antibody Pair -  
BSA and Azide free (ab270374)

Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab269372](#)), which uses the same antibody pair. For additional information on the performance of pair and kit, refer to the corresponding kit datasheet. Due to differences in the formulation and format of the antibodies in this pair, they cannot be used as substitutes for the antibody components in our SimpleStep ELISA® Kits.

Powered by  
recombinant antibodies



**Research with confidence**  
Consistent and reproducible results



**Long-term and scalable supply**  
Recombinant technology



**Success from the first experiment**  
Confirmed specificity



**Ethical standards compliant**  
Animal-free production

Sandwich ELISA - Human Fetuin-A Antibody Pair -  
BSA and Azide free (ab270374)

To learn more about the advantages of recombinant antibodies see [here](#).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors