abcam

Product datasheet

Anti-HSV1 + HSV2 antibody ab226858

2 References

製品の概要

製品名 Anti-HSV1 + HSV2 antibody

製品の詳細 Rabbit polyclonal to HSV1 + HSV2

由来種 Rabbit

特異性 ab226858 is specific for HSV1 and HSV2. It recognizes ICPs and late structural (virion) antigens.

アプリケーション **適用あり**: WB

種交差性 交差種: Herpes simplex virus, Human herpesvirus 2

免疫原 Tissue, cells or virus corresponding to Herpes simplex virus HSV1 + HSV2. Strain F.

特記事項 Concentration is typically 4-5 mg/ml.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years.

Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets

your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be

found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long

term. Avoid freeze / thaw cycle.

バッファー pH: 7.20

Preservative: 0.1% Sodium azide

Constituent: PBS

特記事項(精製) Purity >95% **ポリ/モノ** ポリクローナル

アイソタイプ lqG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab226858の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration.

ターゲット情報

関連性

There are two main types of Herpes Simplex Virus (HSV), 1 and 2. The primary difference between the two types is in where they typically establish latency in the body. HSV1 usually establishes latency in the trigeminal ganglion, a collection of nerve cells near the ear. From there, it tends to recur on the lower lip or face causing cold sores. HSV2 usually resides in the sacral ganglion at the base of the spine. From there, it recurs in the genital area but even this difference is not absolute as in rare cases it can infect any part of the body, including the eyes. Up to 30% of genital herpes are also caused by HSV1. Both HSV1 and 2 are highly contagious through contact.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors