# abcam

## Product datasheet

## Recombinant Human Dopamine D2 Receptor protein ab112281

## 画像数1

製品の詳細

製品名 Recombinant Human Dopamine D2 Receptor protein

生理活性 useful for Antibody Production and Protein Array

**発現系** Wheat germ

アクセッション番号 <u>P14416</u>

タンパク質長 Protein fragment

Animal free No.

由来 Recombinant

生物種 Human

配列 MDPLNLSWYDDDLERQNWSRPFNGSDGKADRPHYNYYATLLT

LLIAVIVF

GNVLVCMAVSREKALQTTTNYLIVSLAVADLLVATLVMPWVV

YLEVVG EW KFSRIHCDIF

予測される分子量 38 kDa including tags

**領域** 1 to 110

タブ GST tag N-Terminus

## 特性

Our Abpromise guarantee covers the use of ab112281 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Western blot

SDS-PAGE

**ELISA** 

Peptide Array

製品の状態 Liquic

備考 This product is useful for Antibody Production and Protein Array.

## 前処理および保存

#### 保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCI

Reduced glutathione

## 関連情報

機能 This is one of the five types (D1 to D5) of receptors for dopamine. The activity of this receptor is

mediated by G proteins which inhibit adenylyl cyclase.

**関連疾患** Defects in DRD2 are associated with dystonia type 11 (DYT11) [MIM:159900]; also known as


alcohol-responsive dystonia. DYT11 is a myoclonic dystonia. Dystonia is defined by the presence of sustained involuntary muscle contractions, often leading to abnormal postures. DYT11 is characterized by involuntary lightning jerks and dystonic movements and postures alleviated by alcohol. Inheritance is autosomal dominant. The age of onset, pattern of body involvement,

presence of myoclonus and response to alcohol are all variable.

配列類似性 Belongs to the G-protein coupled receptor 1 family.

細胞内局在 Cell membrane.

#### 画像



12.5% SDS-PAGE Stained with Coomassie Blue with recombinant protein.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <a href="https://www.abcam.co.jp/abpromise">https://www.abcam.co.jp/abpromise</a> or contact our technical team.

## Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors