abcam

Product datasheet

20nm Gold Nanoparticles (10 OD) ab269936

画像数3

製品の概要

製品名

20nm Gold Nanoparticles (10 OD)

特記事項

20nm Gold Nanoparticles (10 OD) (ab269936) is an aqueous suspension of spherical metallic nanoparticles manufactured to the highest standard to ensure:

Uniform spherical shape

Narrow size distribution: 20nm ± 2nm, <10% CV

Ultra high quality

Batch to batch consistency

High efficiency antibody binding for maximum sensitivity

Available at high concentrations to simplify assay development.

This product is manufactured by Expedeon, an Abcam company. It was previously called 20nm Gold Nanoparticles (10 OD). Expedeon product code 3210-0010 is the same as the 10 mL size of this product and code 3210-0100 is the same as the 100 mL size.

Our colloidal Gold Nanoparticles are suspended in citrate buffer and are suitable for:

Protein adsorption

Oligonucleotide adsorption Surface modifications.

Narrow size distribution and consistent shape is important and preferable for most applications.

A consistent, spherical shape means that particles are less prone to aggregation, which can be a problem associated with gold nanoparticles. Furthermore uniform shape ensures consistent binding of antibodies, ensuring reliable assay results.

製品の特性

製品の状態 Liqu


保存方法 Shipped at Room Temperature. Store at +4°C. Do Not Freeze.

パッファー pH: 7.00

Constituents: 0.0566% Gold nanoparticles, 0.00003% Sodium citrate


画像

1





Electron Microscopy - 20nm Gold Nanoparticles (10 OD) (ab269936)

Electron microscopy (TEM) example pictures showing the high consistency of size and shape of 20nm Gold Nanoparticles (10 OD) (ab269936) (A) compared to that of aggregated, irregular shaped and sized gold nanoparticles of a leading competitor (B).



Spectrum - 20nm Gold Nanoparticles (10 OD) (ab269936)





Functional Studies - 20nm Gold Nanoparticles (10 OD) (ab269936)

Size distribution of 20nm Gold Nanoparticles (10 OD) (ab269936) displaying a narrow size distribution with a CV (standard deviation / diameter) of less than 10%.

Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

· Guarantee only valid for products bought direct from Abcam or one of our authorized distributors