abcam

Product datasheet

20x TBS-T with Tween 20 ab64204

★★★★★ 1 Abreviews 1 References 画像数 1

製品の概要

製品名 20x TBS-T with Tween 20

アプリケーション 適用あり: IHC-P

特記事項 Tris Buffered Saline (TBS) wash buffer with Tween 20 for use during immunolabeling techniques

such as IHC.

The salts in 20x solution are at a very high concentration. It has been found that at temperatures below 44 degrees F (6 degrees C), the salts fall out of solution. If this occurs during your shipment

conditions, please stir the product to get the salts back into solution.

製品の特性

製品の状態 Liquid

保存方法 Store at room temperature.

バッファー pH: 7.2

Constituents: 99% Tris buffered saline (x20), 1% Tween

20x Concentrated

関連性 TBS is a solution of Tris Buffered Saline for use with immunolabeling techniques as a wash buffer.

Tween 20 is sometimes added to promote more effective washings, resulting in decreased non-

specific background staining. Tween 20 is also known to promote reagent spreading in


automated staining methodologies.

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab64204の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent dilution.



Immunohistochemistry (Formalin/PFA-fixed paraffinembedded sections) - 20x TBS-T with Tween 20 (ab64204)

ab64204 20X TBS with Tween 20 used along with an anti-CD20 antibody on Human Tonsil

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- · Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit https://www.abcam.co.jp/abpromise or contact our technical team.

Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors