

Extraction Buffer 5X PTR ab193970

SimpleStep ELISA[®]

16 References

製品の概要

製品名 Extraction Buffer 5X PTR

製品の概要 Abcam's Extraction Buffer 5X is formulated for the preparation of cell or tissues extracts for use in immunoassays.

The Extraction Buffer should be diluted to 1X in the appropriate buffer before use. If required the Extraction Buffer 5X can be used in conjunction of Abcam's Extraction Enhancer Buffer 50X (**ab193971**).

特記事項

Instructions for use

Before use the Extraction Buffer 5X should be diluted to 1X with distilled water.

Note: The provided Extraction Buffer 5X contains phosphatase inhibitors and protease inhibitor aprotinin. Additional protease inhibitors can be added if required.

If the Extraction Buffer 5X is to be used in conjunction with Extraction Enhancer Buffer 50X (**ab193971**), then follow the instructions below:

1X Extraction Buffer + Enhancer

Prepare 1X Extraction Buffer + Enhancer by diluting Extraction Buffer 5X and Extraction Enhancer Buffer 50X to 1X with deionized water. To make 10 mL 1X Extraction Buffer + Enhancer combine 7.8 mL deionized water, 2 mL Extraction Buffer 5X and 200 µL Extraction Enhancer Buffer 50X. Mix thoroughly and gently.

Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances.

It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.

製品の特性

保存方法 Store at +4°C. Please refer to protocols.

内容	10 ml	50 ml
5X Cell Extraction Buffer PTR (ab193970)	1 x 10ml	1 x 50ml

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors