

## Product datasheet

# 125ml 25x TBS (pH 7.4) ab64203

### 製品の概要

製品名	125ml 25x TBS (pH 7.4)
アプリケーション	適用あり: IHC-P
特記事項	In order to provide this product at 25X, the salts in solution are at a very high concentration. It has been found that at temperatures below 44°F (6°C), the salts fall out of solution. If this occurs during your shipment conditions, please stir the product to get the salts back into solution.

### 製品の特性

製品の状態	Liquid
保存方法	Store at room temperature.
バッファー	pH: 7.40 Constituent: 100% Tris buffered saline (x25) Note: Buffer 25x concentrated.
関連性	TBS is a solution of Tris Buffered Saline for use with immunolabeling techniques as a wash buffer. Tween 20 is sometimes added to promote more effective washings, resulting in decreased non-specific background staining. Tween 20 is also known to promote reagent spreading in automated staining methodologies.

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab64203** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors