

Product datasheet

125ml 25x TBS (pH 7.4) ab64203

製品の概要

製品名	125ml 25x TBS (pH 7.4)
アプリケーション	適用あり: IHC-P
特記事項	In order to provide this product at 25X, the salts in solution are at a very high concentration. It has been found that at temperatures below 44°F (6°C), the salts fall out of solution. If this occurs during your shipment conditions, please stir the product to get the salts back into solution.

製品の特性

製品の状態	Liquid
保存方法	Store at room temperature.
バッファー	pH: 7.40 Constituent: 100% Tris buffered saline (x25) Note: Buffer 25x concentrated.
関連性	TBS is a solution of Tris Buffered Saline for use with immunolabeling techniques as a wash buffer. Tween 20 is sometimes added to promote more effective washings, resulting in decreased non-specific background staining. Tween 20 is also known to promote reagent spreading in automated staining methodologies.

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab64203** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors