

Product datasheet

Anti-6X His tag® antibody [HIS-1] ab49936

製品の概要

製品名	Anti-6X His tag® antibody [HIS-1]
製品の詳細	Mouse monoclonal [HIS-1] to 6X His tag®
由来種	Mouse
アプリケーション	適用あり: Flow Cyt, WB, Dot blot, ELISA, ICC/IF
免疫原	Recombinant full length protein corresponding to 6X His tag®.
特記事項	<p>This product was changed from ascites to tissue culture supernatant on 25th October 2016. The following lot(s) is/are from ascites and is still in stock as of 25th October 2016 - GR224242. Lot numbers higher than GR224242 will be from tissue culture supernatant. Please note that the dilutions may need to be adjusted accordingly.</p> <p>HIS-TAG® is a trademark of EMD Biosciences, Inc.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: PBS</p>
精製度	Proprietary Purification
特記事項(精製)	Purified from Tissue culture supernatant.
ポリ/モノ	モノクローナル
クローン名	HIS-1
アイソタイプ	IgG2a

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab49936** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration.
WB		1/3000.
Dot blot		Use at an assay dependent concentration. Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.

ターゲット情報

関連性 The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.

細胞内局在 Depends upon the localization of the parent protein tagged with hexahistidine.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors