

Product datasheet

Anti-6X His tag® antibody (HRP) ab3553

★★★★★ 1 Abreviews 1 References

製品の概要

製品名	Anti-6X His tag® antibody (HRP)
製品の詳細	Chicken polyclonal to 6X His tag® (HRP)
由来種	Chicken
標識	HRP
アプリケーション	適用あり: WB, ELISA
免疫原	Synthetic peptide: HHHHHH conjugated to KLH.
特記事項	HIS-TAG® is a trademark of EMD Biosciences, Inc.

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	50 mM HEPES pH 7.0, 0.1 M NaCl, 0.2% BSA, 0.01% Thimerosal
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
ミエローマ	unknown
アイソタイプ	IgY
軽鎖の種類	unknown

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab3553** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリ ケーショ ン	Abreviews	特記事項
WB	★★★★★	Use at an assay dependent concentration. 1/1000.
ELISA		Use at an assay dependent concentration. a 1/10000 dilution gave an O.D.= 1.0 in a 15 minute reaction against peptide conjugated to a different carrier than used for anti-peptide purification.

ターゲット情報

関連性	The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. The antibody recognizes the His-tag fused to the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.
細胞内局在	Depends upon the localization of the parent protein tagged with hexahistidine.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors