

Product datasheet

Anti-ABCA6 antibody ab100859

画像数 1

製品の概要

製品名	Anti-ABCA6 antibody
製品の詳細	Rabbit polyclonal to ABCA6
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	16 amino acids synthetic peptide derived from the N terminal region of Human ABCA6.
ポジティブ・コントロール	Fetal brain lysate; 293T cell lysate.

法規制情報

医薬用外毒物

製品の特性

製品の状態	Lyophilised:Reconstitute in 200ul sterile water
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituent: 1% BSA
精製度	Affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab100859** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 184 kDa (predicted molecular weight: 184 kDa).

ターゲット情報

機能	Probable transporter which may play a role in macrophage lipid homeostasis.
組織特異性	Widely expressed with higher expression in liver.
配列類似性	Belongs to the ABC transporter superfamily. ABCA family. Contains 2 ABC transporter domains.
発生段階	Expressed in fetal kidney, lung and liver.
細胞内局在	Membrane.

画像



All lanes : Anti-ABCA6 antibody (ab100859)
at 1/1000 dilution

Lane 1 : Human fetal brain lysate

Lane 2 : 293T cell lysate

Predicted band size: 184 kDa

Observed band size: 184 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors