

Product datasheet

Anti-ANXA9 antibody ab37501

画像数 2

製品の概要

製品名	Anti-ANXA9 antibody
製品の詳細	Chicken polyclonal to ANXA9
由来種	Chicken
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human
免疫原	Full length native protein (purified), corresponding to amino acids 1-338 of Human ANXA9

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab37501** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use a concentration of 1 µg/ml.
WB		1/500. Predicted molecular weight: 38 kDa.

ターゲット情報


機能	Low affinity receptor for acetylcholine known to be targeted by disease-causing pemphigus vulgaris antibodies in keratinocytes.
組織特異性	Expressed in the stratified squamous skin epithelium, but not in epithelia of other types (at protein level).
配列類似性	Belongs to the annexin family. Contains 4 annexin repeats.

画像



ELISA - Anti-ANXA9 antibody (ab37501)

Direct ELISA test. Native Human Annexin A9 was test antigen. ab37501 (chicken IgY anti-annexin A9) was primary antibody and Goat anti-IgY-HRP was secondary antibody. A fixed amount of antibody (1 µg/ml) and serial dilutions of antigen were used.



Western blot - Anti-ANXA9 antibody (ab37501)

Anti-ANXA9 antibody (ab37501) at 1/2000 dilution + E. coli-derived fusion protein

Secondary

Goat anti-IgY-HRP at 1/1000 dilution

Predicted band size: 38 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors