

## Product datasheet

# Anti-AP3 antibody ab85075

### 画像数 1

#### 製品の概要

---

製品名	Anti-AP3 antibody
製品の詳細	Rabbit polyclonal to AP3
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	A 15 amino acid peptide from near the C terminal residues of Human AP3, swissprot O14617
ポジティブ・コントロール	HeLa cell lysate.

#### 法規制情報

---

### 医薬用外毒物

#### 製品の特性

---

製品の状態	Lyophilised:Reconstitute with 200ul distilled sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium Azide
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab85075** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/200 - 1/1000. Detects a band of approximately 130 kDa (predicted molecular weight: 130 kDa).

## ターゲット情報

機能	Part of the AP-3 complex, an adaptor-related complex which is not clathrin-associated. The complex is associated with the Golgi region as well as more peripheral structures. It facilitates the budding of vesicles from the Golgi membrane and may be directly involved in trafficking to lysosomes.
組織特異性	Present in all adult tissues examined with the highest levels in skeletal muscle, heart, pancreas and testis.
配列類似性	Belongs to the adaptor complexes large subunit family. Contains 11 HEAT repeats.
細胞内局在	Cytoplasm. Golgi apparatus membrane.

## 画像



Anti-AP3 antibody (ab85075) at 1/500 dilution  
+ HeLa cell lysate

**Predicted band size:** 130 kDa

**Observed band size:** 130 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors