

Product datasheet

Anti-APRIL antibody (Biotin) ab84029

製品の概要

製品名	Anti-APRIL antibody (Biotin)
製品の詳細	Rabbit polyclonal to APRIL (Biotin)
標識	Biotin
アプリケーション	適用あり: WB, ELISA, Sandwich ELISA
種交差性	交差種: Human
免疫原	Highly pure (>98%) recombinant full length APRIL protein (Human).

製品の特性

製品の状態	Lyophilised:Centrifuge vial prior to opening. Reconstitute in 50 µl sterile PBS containing 0.1% BSA to give a concentration of 1.0 mg/ml.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: 0.1% BSA, PBS, pH 7.2
精製度	Protein A purified
特記事項(精製)	ab84029 is a sterile filtered antibody solution.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab84029** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.1 - 0.2 µg/ml. Predicted molecular weight: 27 kDa.
ELISA		Use a concentration of 0.25 - 1 µg/ml.
Sandwich ELISA		Use a concentration of 0.25 - 1 µg/ml.

ターゲット情報

機能	Cytokine that binds to TNFRSF13B/TACI and to TNFRSF17/BCMA. May be implicated in the regulation of tumor cell growth. May be involved in monocyte/macrophage-mediated immunological processes.
組織特異性	Expressed at high levels in transformed cell lines, cancers of colon, thyroid, lymphoid tissues and specifically expressed in monocytes and macrophages.
配列類似性	Belongs to the tumor necrosis factor family.
翻訳後修飾	The precursor is cleaved by furin.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors