

Product datasheet

Anti-ARHGEF7 antibody ab125870

画像数 2

製品の概要

製品名	Anti-ARHGEF7 antibody
製品の詳細	Rabbit polyclonal to ARHGEF7
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to a region within amino acids 218-428 of Human ARHGEF7 (Q14155).
ポジティブ・コントロール	NT2D1, IMR32, U-87 MG and MCF7 whole cell lysates.
特記事項	Keep as concentrated solution

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 78% PBS, 20% Glycerol, 1% BSA
精製度	Immunogen affinity purified
特記事項(精製)	ab125870 is purified by antigen-affinity chromatography.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab125870** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 90 kDa.
IHC-P		1/100 - 1/1000. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. Use 10mM Citrate buffer (pH6.0) or Tris-EDTA buffer (pH8.0).

ターゲット情報

機能	Acts as a RAC1 guanine nucleotide exchange factor (GEF) and can induce membrane ruffling. May function as a positive regulator of apoptosis. May function in cell migration.
配列類似性	Contains 1 CH (calponin-homology) domain. Contains 1 DH (DBL-homology) domain. Contains 1 PH domain. Contains 1 SH3 domain.


画像



Anti-ARHGEF7 antibody (ab125870) at
1/1000 dilution + IMR32 whole cell lysate at 30
µg

Predicted band size: 90 kDa

7.5% SDS PAGE



ab125870, at 1/500, staining ARHGEF7 in paraffin embedded Human AdCA (adenocarcinoma) by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ARHGEF7 antibody (ab125870)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors