

Product datasheet

Anti-Alpha Taxilin antibody ab37882

画像数 1

製品の概要

製品名	Anti-Alpha Taxilin antibody
製品の詳細	Chicken polyclonal to Alpha Taxilin
由来種	Chicken
アプリケーション	適用あり: WB
種交差性	交差種: Dog, Human, Chimpanzee
免疫原	Synthetic peptide, corresponding to amino acids 409-539 of Human Alpha Taxilin.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab37882** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/2000. Detects a band of approximately 65 kDa (predicted molecular weight: 62 kDa).

ターゲット情報

機能	May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in
----	---

neuroendocrine cells.


組織特異性

Ubiquitous, with much higher expression in heart, kidney, liver and pancreas.

配列類似性

Belongs to the taxilin family.

画像



Anti-Alpha Taxilin antibody (ab37882) at
1/2000 dilution + E. coli derived fusion protein
at 2.5 µg

Secondary

Goat anti IgY HRP at 1/1000 dilution

Predicted band size: 62 kDa

Observed band size: 65 kDa

The arrow next to the image points to the
Alpha Taxilin fusion protein.

Colorimetric method was used for signal
development.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors