

## Product datasheet

# Anti-BRD1 antibody [2093C5a] ab50989

### 画像数 1

#### 製品の概要

製品名	Anti-BRD1 antibody [2093C5a]
製品の詳細	Mouse monoclonal [2093C5a] to BRD1
由来種	Mouse
アプリケーション	<b>適用あり:</b> Dot blot, IP, WB
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment (Human)
ポジティブ・コントロール	RIPA lysate of HeLa cells.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
精製度	Protein G purified
特記事項(精製)	Filtered through a 0.22 µm membrane.
ポリモノ	モノクローナル
クローン名	2093C5a
アイソタイプ	IgG1

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab50989** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
IP		Use at 10 µg/mg of lysate.
WB		1/50. Predicted molecular weight: 137 kDa.

## ターゲット情報

機能	Component of the MOZ/MORF complex which has a histone H3 acetyltransferase activity.
組織特異性	Highly expressed in testis.
配列類似性	Contains 1 bromo domain. Contains 1 PHD-type zinc finger. Contains 1 PWWP domain.
細胞内局在	Nucleus.

## 画像



Immunoprecipitation analysis of RIPA lysate of HeLa cells at 500 µg using ab50989 (1/50).  
Secondary detection using anti mouse IgG (1/2500).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors