

## Product datasheet

# Anti-BXDC2 antibody [2092C1a] ab51271

### 画像数 1

#### 製品の概要

製品名	Anti-BXDC2 antibody [2092C1a]
製品の詳細	Mouse monoclonal [2092C1a] to BXDC2
由来種	Mouse
アプリケーション	<b>適用あり:</b> Dot blot, WB
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment (Human)

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
精製度	Protein G purified
特記事項 (精製)	ab51271 was purified from culture supernatant of the hybridoma, cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum.
ポリ/モノ	モノクローナル
クローン名	2092C1a
アイソタイプ	IgG1

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab51271** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 41 kDa.

## ターゲット情報


---

機能	Required for biogenesis of the 60S ribosomal subunit.
配列類似性	Belongs to the BRX1 family. Contains 1 Brix domain.
細胞内局在	Nucleus > nucleolus.

---

## 画像

---



Anti-BXDC2 antibody [2092C1a] (ab51271)

**Predicted band size:** 41 kDa

**Observed band size:** 38 kDa

Western blot analysis of immunised recombinant protein  
Western blot analysis of immunised recombinant protein

Western blot - Anti-BXDC2 antibody [2092C1a] (ab51271)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors