

Product datasheet

Anti-Bex1 antibody ab107215

画像数 1

製品の概要

製品名	Anti-Bex1 antibody
製品の詳細	Rabbit polyclonal to Bex1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to internal sequence amino acids 66-94 of Human Bex1 (NP_060946.3)
ポジティブ・コントロール	WiDr cell line lysates

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	This antibody is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab107215** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 15 kDa.

ターゲット情報

関連性	Bex1 is a signaling adapter molecule involved in p75NTR/NGFR signaling. It is thought to play a role in cell cycle progression and neuronal differentiation. Bex1 inhibits neuronal differentiation in response to nerve growth factor (NGF). Bex1 has also been implicated as a tumor suppressor
細胞内局在	Cytoplasmic and Nuclear

画像



Anti-Bex1 antibody (ab107215) at 1/100 dilution + WiDr cell line lysates at 35 µg

Predicted band size: 15 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors