

Product datasheet

Anti-Blood Group A1B Antigen antibody [HE-24] ab2525

2 References

製品の概要

製品名	Anti-Blood Group A1B Antigen antibody [HE-24]
製品の詳細	Mouse monoclonal [HE-24] to Blood Group A1B Antigen
由来種	Mouse
特異性	This antibody recognises human blood group A ₁ B antigen. It does not detect A ₂ B antigen.
アプリケーション	適用あり: Agglutination
種交差性	交差種: Human
免疫原	Glycoprotein fraction isolated from the saliva of secretors with blood group A.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: Tissue Culture Supernatant
精製度	Tissue culture supernatant
特記事項 (精製)	Concentrated by ultrafiltration (100 kDa cut-off). Actual immunoglobulin concentration not determined.
ポリ/モノ	モノクローナル
クローン名	HE-24
ミエローマ	unknown
アイソタイプ	IgM
軽鎖の種類	unknown

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab2525** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Agglutination		Use at an assay dependent concentration. The minimum titre is 1:128 with erythrocytes of blood group A ₁ B.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors