

## Product datasheet

# Anti-C2CD2 antibody ab101445

### 画像数 1

#### 製品の概要

---

製品名	Anti-C2CD2 antibody
製品の詳細	Rabbit polyclonal to C2CD2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to a region within the amino acids 1-165 of Human C2CD2 (AAH62323).
ポジティブ・コントロール	NT2D1, IMR32 and MCF-7.

#### 法規制情報

---

### 医薬用外毒物

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab101445** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

## アプリケーション

## Abreviews

## 特記事項

WB


1/500 - 1/3000. Predicted molecular weight: 42 kDa.

## ターゲット情報

細胞内局在

Secreted

## 画像



Anti-C2CD2 antibody (ab101445) at 1/1000  
dilution + IMR32 whole cell lysate at 30 µg

**Predicted band size:** 42 kDa

10% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors