

Product datasheet

Anti-C5a / C5a des Arg antibody [2952] ab11878

★★★★☆ 1 Abreviews 3 References 画像数 1

製品の概要

製品名	Anti-C5a / C5a des Arg antibody [2952]
製品の詳細	Mouse monoclonal [2952] to C5a / C5a des Arg
特異性	This antibody reacts with a neo-epitope (des-Arg) on C5a that is formed upon cleavage of C5 into C5a and C5b. It will detect C5a and C5a des-Arg (Not C5).
アプリケーション	適用あり: IHC-Fr, Functional Studies, IHC-P, ELISA, WB
種交差性	交差種: Human
免疫原	ISHKDMQLG (C5a amino acid 65-73)
	Run BLAST with Run BLAST with
ポジティブ・コントロール	Human Granulocytes

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 0.1% BSA, PBS
精製度	Protein G purified
特記事項(精製)	0.2 µm filtered
ポリ/モノ	モノクローナル
クローン名	2952
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab11878** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.
Functional Studies		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
ELISA	★★★★☆	Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

ターゲット情報

関連性 C5 is involved in the activation of the lytic pathway within the complement system, which is an important factor in innate immunity.

画像



ELISA with [ab11876](#) as capture antibody at 1/100, [ab11878](#) as detection antibody at 1/250, diluted in carbonate buffer or PBS. Dilutions of commercial C5a were used as the standard.

ELISA - Anti-C5a / C5a des Arg antibody [2952]
([ab11878](#))

This image is courtesy of an anonymous Abreview

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors