

Product datasheet

Anti-C9orf98 antibody ab95948

画像数 2

製品の概要

製品名	Anti-C9orf98 antibody
製品の詳細	Rabbit polyclonal to C9orf98
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, ICC/IF
種交差性	<b>交差種:</b> Human
免疫原	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 242-461 of C9orf98 (Human)
ポジティブ・コントロール	293T; A431; H1299; HeLa S3; HepG2; Molt 4 or Raji cell lines; A431 and HeLa S3 whole cell lysate (WB); HeLa cells (IF)

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab95948** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 55 kDa.
ICC/IF		1/100 - 1/200.

## ターゲット情報

**関連性** Belongs to the adenylate kinase family.

## 画像




**All lanes** : Anti-C9orf98 antibody (ab95948)  
at 1/500 dilution

**Lane 1** : A431 whole cell lysate  
**Lane 2** : HeLa S3 whole cell lysate

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 55 kDa

gel concentration: 7.5%



Immunofluorescence analysis of  
paraformaldehyde fixed HeLa cells, using  
ab95948 at 1/200 dilution:

- (1) ab95948
- (2) ab95948 merged with DNA probe.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors