

## Product datasheet

### Anti-CCR1 antibody [MM0061-7B17] ab89055

#### 製品の概要

---

製品名	Anti-CCR1 antibody [MM0061-7B17]
製品の詳細	Mouse monoclonal [MM0061-7B17] to CCR1
由来種	Mouse
特異性	ab89055 detects CCR1, and not other CCRs.
アプリケーション	<b>適用あり:</b> IHC-P, Flow Cyt
種交差性	<b>交差種:</b> Human
免疫原	293 cells transfected with human CCR1.
特記事項	This antibody was produced from a hybridoma (mouse myeloma fused with spleen cells from an immunized mouse).

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS
精製度	Protein G purified
特記事項 (精製)	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and lyophilized from a 0.2 µm filtered solution.
ポリ/モノ	モノクローナル
クローン名	MM0061-7B17
アイソタイプ	IgG2

#### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab89055** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

IHC-P

Flow Cyt

**追加情報**

Flow Cyt: 1/100 - 1/800.  
IHC-P: 1/50 - 1/200.

Not yet tested in other applications.  
Optimal dilutions/concentrations should be determined by the end user.

**ターゲット情報**

**機能**

Receptor for a C-C type chemokine. Binds to MIP-1-alpha, MIP-1-delta, RANTES, and MCP-3 and, less efficiently, to MIP-1-beta or MCP-1 and subsequently transduces a signal by increasing the intracellular calcium ions level. Responsible for affecting stem cell proliferation.

**組織特異性**

Widely expressed in different hematopoietic cells.

**配列類似性**

Belongs to the G-protein coupled receptor 1 family.

**細胞内局在**

Cell membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors