

Product datasheet

Anti-CD15 antibody [28] ab20137

3 References

製品の概要

| | |
|----------|--|
| 製品名 | Anti-CD15 antibody [28] |
| 製品の詳細 | Mouse monoclonal [28] to CD15 |
| 特異性 | Specific to CD15. 28 is useful for leukaemia typing and purification of stem cells. |
| アプリケーション | 適用あり: IP, Flow Cyt, IHC-Fr |
| 種交差性 | 交差種: Human |
| 免疫原 | Cell preparation: Monocytes, separated from other peripheral blood leucocytes (PBL) on fibronectin plates. |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: None Constituents: PBS |
| 精製度 | IgG fraction |
| ポリ/モノ | モノクローナル |
| クローン名 | 28 |
| ミエローマ | Sp2/0-Ag14 |
| アイソタイプ | IgM |
| 軽鎖の種類 | unknown |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab20137** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| IP | | |
| Flow Cyt | | |

IHC-Fr

追加情報

Flow Cyt: Use at an assay dependent dilution.

IHC-Fr: Use at an assay dependent dilution.

Note: samples should be acetone fixed.

IP: Use at an assay dependent dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報**関連性**

CD15 is a carbohydrate adhesion molecule (and not a protein) that mediates phagocytosis and chemotaxis. Synthesis is directed by FUT4 in lymphoid cells and mature granulocytes, and by FUT9 in promyelocytes and monocytes.

細胞内局在

Golgi Apparatus; Membrane-bound form in trans cisternae of Golgi.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors