

Product datasheet

Anti-CD23 antibody [B-G6] (Phycoerythrin) ab27304

画像数 1

製品の概要

製品名	Anti-CD23 antibody [B-G6] (Phycoerythrin)
製品の詳細	Mouse monoclonal [B-G6] to CD23 (Phycoerythrin)
由来種	Mouse
標識	Phycoerythrin. Ex: 488nm, Em: 575nm
特異性	This antibody is specific for CD23
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Daudi cell line

法規制情報

医薬用外毒物

製品の特性

製品の状態	Lyophilised:Reconstitute with 1 ml deionised water
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, Phosphate-buffered saline
精製度	Ion Exchange Chromatography
ポリ/モノ	モノクローナル
クローン名	B-G6
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab27304** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Flow Cyt

追加情報

Flow Cyt: Use 10 μ l for 10⁶ cells or 100 μ l of whole blood.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能

This receptor has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).

配列類似性

Contains 1 C-type lectin domain.


翻訳後修飾

N- and O-glycosylated.

細胞内局在

Cell membrane. Secreted. Also exists as a soluble excreted form.

画像



A typical staining pattern with ab27304 (10 μ l per 10⁶ cells or 100 μ l of whole blood) of lymphocytes.

Flow Cytometry - Anti-CD23 antibody [B-G6]
(Phycoerythrin) (ab27304)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors