

Product datasheet

Anti-CD23 antibody [EBVCS2] (Allophycocyanin) ab95573

★★★★☆ 1 Abreviews [画像数 1](#)

製品の概要

製品名	Anti-CD23 antibody [EBVCS2] (Allophycocyanin)
製品の詳細	Mouse monoclonal [EBVCS2] to CD23 (Allophycocyanin)
由来種	Mouse
標識	Allophycocyanin. Ex: 645nm, Em: 660nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Epstein-Barr Virus-transformed Human cells
ポジティブ・コントロール	Normal Human peripheral blood cells

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS, 150mM Sodium chloride, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	EBVCS2
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab95573** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt	★★★★☆	Use 5µl for 10 ⁶ cells.

ターゲット情報

機能	This receptor has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).
配列類似性	Contains 1 C-type lectin domain.
翻訳後修飾	N- and O-glycosylated.
細胞内局在	Cell membrane. Secreted. Also exists as a soluble excreted form.

画像



Flow Cytometry - Anti-CD23 antibody [EBVCS2]
(Allophycocyanin) (ab95573)

Staining of normal human peripheral blood cells with Anti-Human CD19 FITC and Mouse IgG1 ? Isotype Control APC (left) or ab95573 (right). Cells in the lymphocyte gate were used for analysis.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors