

Product datasheet

Anti-CD40 antibody [HI40a] (CF405M) ab119494

画像数 1

製品の概要

製品名	Anti-CD40 antibody [HI40a] (CF405M)
製品の詳細	Mouse monoclonal [HI40a] to CD40 (CF405M)
由来種	Mouse
標識	CF405M. Ex: 408nm, Em: 452nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Tonsil cells.
ポジティブ・コントロール	Normal Human peripheral blood cells.
特記事項	CF405M (Abs/Em Max: 408/450nm). Direct replacement for Pacific Blue dye [®] , BD Horizon™ V450.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Note: Buffer containing antibody stabilizer solution.
精製度	Immunogen affinity purified
ポリ/モノ	モノクローナル
クローン名	HI40a
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab119494** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors