

## Product datasheet

# Anti-CD99 antibody [HI156] (Phycoerythrin) ab91298

### 画像数 2

#### 製品の概要

製品名	Anti-CD99 antibody [HI156] (Phycoerythrin)
製品の詳細	Mouse monoclonal [HI156] to CD99 (Phycoerythrin)
標識	Phycoerythrin. Ex: 488nm, Em: 575nm
特異性	ab91298 recognizes the MIC 2 (E2) antigen expressed on most cells except those of the myeloid family.
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Human
免疫原	Purified E rosette forming cells from Human PBL.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 1% BSA, pH 7.2
精製度	Immunogen affinity purified
ポリ/モノ	モノクローナル
クローン名	HI156
アイソタイプ	IgG2a

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab91298** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

#### 追加情報

Flow Cyt: Use 20 µl for 10<sup>6</sup> cells.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


機能	Involved in T-cell adhesion processes. It is involved in spontaneous rosette formation with erythrocytes.
配列類似性	Belongs to the CD99 family.
翻訳後修飾	Extensively O-glycosylated.
細胞内局在	Membrane.

## 画像



ab91298, at 20  $\mu\text{l}/10^6$  cells, staining CD99 in normal Human peripheral blood by Flow cytometry.

Flow Cytometry - CD99 antibody [HI156]  
(Phycoerythrin) (ab91298)



ab91298, at 20  $\mu\text{l}/10^6$  cells, staining CD99 in normal Human peripheral blood by Flow cytometry.

Flow Cytometry - CD99 antibody [HI156]  
(Phycoerythrin) (ab91298)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors