

Product datasheet

Anti-CTAG1B antibody ab113918

画像数 2

製品の概要

製品名	Anti-CTAG1B antibody
製品の詳細	Rabbit polyclonal to CTAG1B
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide, corresponding to a region within C terminal amino acids 131-180 (NILTRLTAA DHRQLQLSIS SCLQQLSLLM WITQCFLPVF LAQPPSGQRR) of Human CTAG1B (NM_001327). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Jurkat cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituents: 98% PBS, 2% Sucrose
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab113918** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB	Use a concentration of 1 µg/ml. Predicted molecular weight: 18 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	--

ターゲット情報

組織特異性

Expressed in testis and ovary and in a wide variety of cancers. Detected in uterine myometrium. Expressed from 18 weeks until birth in human fetal testis. In the adult testis, is strongly expressed in spermatogonia and in primary spermatocytes, but not in post-meiotic cells or in testicular somatic cells (at protein level).


配列類似性

Belongs to the CTAG family.

細胞内局在

Cytoplasm.


画像



Anti-CTAG1B antibody (ab113918) at 1 μ g/ml
+ Jurkat cell lysate at 10 μ g

Predicted band size: 18 kDa

Gel concentration 12%



Anti-CTAG1B antibody (ab113918) at 1 μ g/ml
+ 786-0 cell lysate

Predicted band size: 18 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors