


Product datasheet

Anti-Calbindin antibody, prediluted ab53642

画像数 1

製品の概要

製品名	Anti-Calbindin antibody, prediluted
製品の詳細	Rabbit polyclonal to Calbindin, prediluted
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken, Hamster, Cow, Xenopus laevis, Orangutan 
免疫原	Synthetic peptide: KMGKKEFNKAFELYD , corresponding to internal sequence amino acids 185-199 of Human Calbindin Run BLAST with Run BLAST with
ポジティブ・コントロール	Rat brain, human cerebellum.

製品の特性

製品の状態	Prediluted
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 1% BSA, 50mM Tris, pH 7.6
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab53642** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		

追加情報

IHC-P: 1/1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at RT for 20 min.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	Buffers cytosolic calcium. May stimulate a membrane Ca(2+)-ATPase and a 3',5'-cyclic nucleotide phosphodiesterase.
配列類似性	Belongs to the calbindin family. Contains 5 EF-hand domains.
ドメイン	This protein has four functional calcium-binding sites; potential sites II and VI have lost affinity for calcium.

Anti-Calbindin antibody, prediluted 画像



Ab53642 at a ready to use dilution staining human cerebellum; formalin fixed, paraffin embedded.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Calbindin antibody (ab53642)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors