

## Product datasheet

# Anti-Calpain 15 antibody ab44204

### 製品の概要

製品名	Anti-Calpain 15 antibody
製品の詳細	Rabbit polyclonal to Calpain 15
由来種	Rabbit
特異性	The antibody binds to Calpain 15, but does not cross react with the other Calpain family members.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	A synthetic peptide based on the carboxy terminal end of the SOL domain of human Calpain 15.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 50% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab44204** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

### 追加情報

WB: 1/1000 when using colorimetric substrates such as BCIP/NBT, and 1/5000 for chemiluminescent substrates. Detects a band of approximately 115 kDa and a series of further cleaved forms. (predicted molecular weight: 117 kDa). Dilution optimised using Chromogenic detection. Not yet tested in other applications. Optimal dilutions/concentrations should be

determined by the end user.

## ターゲット情報

---

組織特異性	Widely expressed with higher expression in brain.
配列類似性	Belongs to the peptidase C2 family. Contains 1 calpain catalytic domain. Contains 5 RanBP2-type zinc fingers.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors