

Product datasheet

Anti-Calpain 15 antibody ab44204

製品の概要

製品名	Anti-Calpain 15 antibody
製品の詳細	Rabbit polyclonal to Calpain 15
由来種	Rabbit
特異性	The antibody binds to Calpain 15, but does not cross react with the other Calpain family members.
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human
免疫原	A synthetic peptide based on the carboxy terminal end of the SOL domain of human Calpain 15.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 50% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab44204** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

追加情報

WB: 1/1000 when using colorimetric substrates such as BCIP/NBT, and 1/5000 for chemiluminescent substrates. Detects a band of approximately 115 kDa and a series of further cleaved forms. (predicted molecular weight: 117 kDa). Dilution optimised using Chromogenic detection. Not yet tested in other applications. Optimal dilutions/concentrations should be

determined by the end user.

ターゲット情報

組織特異性	Widely expressed with higher expression in brain.
配列類似性	Belongs to the peptidase C2 family. Contains 1 calpain catalytic domain. Contains 5 RanBP2-type zinc fingers.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors