

Product datasheet

Anti-Calretinin antibody [EPR1798(I)] ab76028

リコンビナント RabMAb®

画像数 1

製品の概要

製品名	Anti-Calretinin antibody [EPR1798(I)]
製品の詳細	Rabbit monoclonal [EPR1798(I)] to Calretinin
由来種	Rabbit
アプリケーション	適用あり: WB, IP 適用なし: Flow Cyt, ICC or IHC-P
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide within Human Calretinin (N terminal). The exact sequence is proprietary.
ポジティブ・コントロール	Fetal brain, mouse brain and rat brain lysates.
特記事項	

Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMab® patents](#)

This product is a recombinant rabbit monoclonal antibody.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	PBS 49%, Sodium azide 0.01%, Glycerol 50%, BSA 0.05%
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	EPR1798(I)
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab76028** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション**Abreviews****特記事項**

WB

IP

追加情報

IP: 1/50.

WB: 1/1000 - 1/5000. Predicted molecular weight: 31 kDa.


Is unsuitable for Flow Cyt, ICC or IHC-P.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報**機能** Calretinin is a calcium-binding protein which is abundant in auditory neurons.**組織特異性** Brain.**配列類似性** Belongs to the calbindin family.

Contains 6 EF-hand domains.

画像

Western blot - Anti-Calretinin antibody [EPR1798(I)]
(ab76028)

All lanes : Anti-Calretinin antibody

[EPR1798(I)] (ab76028) at 1/5000 dilution

Lane 1 : fetal brain lysate**Lane 2** : mouse brain lysate**Lane 3** : rat brain lysate

Lysates/proteins at 10 µg per lane.

Secondary**All lanes** : goat anti-rabbit HRP at 1/2000 dilution**Predicted band size:** 31 kDa**Observed band size:** 31 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors