

## Product datasheet

# Anti-Cholecystokinin 8 antibody [27F3..0,5D10] ab37274

### 3 References

#### 製品の概要

製品名	Anti-Cholecystokinin 8 antibody [27F3..0,5D10]
製品の詳細	Mouse monoclonal [27F3..0,5D10] to Cholecystokinin 8
特異性	Reacts with CCK-8 and the other CCK analogs with identical C-terminal ends as well as Gastrin 17. No reactivity was observed with Human Heparin cofactor II.
アプリケーション	<b>適用あり:</b> ELISA, IHC-FoFr
種交差性	<b>交差種:</b> Rat
免疫原	CCK-8 fragment 26-33 GGDRDY(SO3H)MGWMDF

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term.
バッファー	Preservative: None Constituents: 0.5M Sodium chloride, 0.01M PBS, pH 7.4
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	27F3..0,5D10
ミエローマ	x63-Ag8.653
アイソタイプ	IgG2a
軽鎖の種類	kappa

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab37274** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項

アプリ ケー ション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration. We recommend coating with CCK-8 and using the antibody in a competitive assay protocol. Increased reactivity with sulphated CCK-8 compared to non-sulphated CCK-8 has been observed.
IHC- FoFr		1/1000. PubMed: 22892424

## ターゲット情報

### 細胞内局在

Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors