

Product datasheet

Choline Acetyltransferase peptide (736-748) ab45678

製品の概要

製品名 Choline Acetyltransferase peptide (736-748)

製品の詳細

由来 Synthetic  
 アミノ酸配列  
 生物種 Human  
 配列 CKEKATRPSQGHQP  
 領域 736 to 748

特性

Our [Abpromise guarantee](#) covers the use of **ab45678** in the following tested applications.  
 The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-Choline Acetyltransferase antibody ([ab27484](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

関連情報

機能 Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.

関連疾患 Myasthenic syndrome, congenital, 6, presynaptic

配列類似性 Belongs to the carnitine/choline acetyltransferase family.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors