

## Product datasheet

# Anti-Chromogranin C antibody [SGII 9G3/3] ab20246

### 4 References

#### 製品の概要

製品名	Anti-Chromogranin C antibody [SGII 9G3/3]
製品の詳細	Mouse monoclonal [SGII 9G3/3] to Chromogranin C
由来種	Mouse
アプリケーション	<b>適用あり:</b> IP, Flow Cyt, IHC-Fr, WB
種交差性	<b>交差種:</b> Mouse, Rat
免疫原	Synthetic peptide: KEENSRENPF (C terminus of Sgll processing fragment p18 Rat)
	<a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: PBS
ポリ/モノ	モノクローナル
クローン名	SGII 9G3/3
ミエローマ	Sp2
アイソタイプ	IgG2b

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab20246** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. <a href="#">ab170192</a> - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.
IHC-Fr		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 71 kDa.

## ターゲット情報

機能	Secretogranin-2 is a neuroendocrine secretory granule protein, which is the precursor for biologically active peptides.
配列類似性	Belongs to the chromogranin/secretogranin protein family.
細胞内局在	Secreted. Neuroendocrine and endocrine secretory granules.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors