

Product datasheet

Cow Troponin T1 peptide ab42256

製品の概要

製品名 Cow Troponin T1 peptide

製品の詳細

由来 Synthetic
 アミノ酸配列
 生物種 Cow
 配列 APPPAEVPEVHEEVH
 領域 16 to 31

特性

Our [Abpromise guarantee](#) covers the use of **ab42256** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking
 製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.
 Preservative: None
 Constituents: 0.1% Acetic acid

関連情報

機能 Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.
 関連疾患 Nemaline myopathy 5
 配列類似性 Belongs to the troponin T family.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors