

Product datasheet

Anti-DDDDK tag antibody [M2] - (Equivalent to FLAG antibodies from Sigma) (DyLight® 488) ab117505

製品の概要

製品名	Anti-DDDDK tag antibody [M2] - (Equivalent to FLAG antibodies from Sigma) (DyLight® 488)
製品の詳細	Mouse monoclonal [M2] to DDDDK tag - (Equivalent to FLAG antibodies from Sigma) (DyLight® 488)
由来種	Mouse
標識	DyLight® 488. Ex: 493nm, Em: 518nm
アプリケーション	適用あり: ICC/IF
免疫原	Synthetic peptide corresponding to DDDDK tag . Sequence: DYKDDDDK
	Run BLAST with Run BLAST with
特記事項	Stable for 1yr.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.20 Constituent: 99% Tris buffered saline
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	M2
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab117505** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		1/500 - 1/5000.

ターゲット情報

関連性 This is a useful tool for the localisation and characterisation of DDDDK tagged proteins.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors